

## Technical Data Sheet

### Mini-Tester PTM100 pH/Chlorine/Bromine

#### Description





PTM100 is a Water-i.d. Pool tester Box for quickly determination of pH, Chlorine and Bromine in pool waters. Included quickly dissolving indicator RAPID tablets are perfectly matching the Pool tester's colour scales for high precision optical comparison.



#### Technical Data

Product Name	Mini Tester PTM100
Manufacturer	Water-i.d. GmbH
Application	Chlorine, Bromine & pH measurement in pool water
Measurement Range	0.5 – 5.0 mg/l mg/l Cl <sub>2</sub> /Br <sub>2</sub> pH 6.8 – 8.0
Reagent Form	Tablets
Packaging Mini Tester	Plastic box (HDPE)
Packaging Tablets	Blister
Stability	10 years of shelf life
Storage Conditions	Cool and dry place
Dimensions (HxWxD)	9.0 x 11.0 x 3.5
Weight [g]	122

#### Features

-  Fast Dissolving Formula
-  High Accuracy
-  Extended Shelf Life
-  Reliable Quality Control

✓ ISO 9001:2015 certified

Water-i.d. GmbH Daimlerstr. 20 D-76344 Eggenstein www.water-id.com info@water-id.com Tel. +49 (0) 721 - 78 20 29 0 Fax +49 (0) 721 78 20 29 11 / 21  
VAT-ID (USt.): DE 271758660 Registered at / Amtsgericht: Mannheim HRB 709390 General manager / Geschäftsführer: Dipl. Betriebswirt Andreas Hock

#### Bankverbindungen / Bank accounts

<b>EUR:</b>	Sparkasse Karlsruhe	<b>IBAN:</b> DE 49 6605 0101 0022607659
<b>EUR:</b>	Volksbank Karlsruhe	<b>IBAN:</b> DE 81 6619 0000 0010 327210
<b>USD only:</b>	Sparkasse Karlsruhe	<b>IBAN:</b> DE 19 6605 0101 0075101204
<b>GBP only:</b>	Sparkasse Karlsruhe	<b>IBAN:</b> DE 18 6605 0101 0075103259
<b>CNY only:</b>	Sparkasse Karlsruhe	<b>IBAN:</b> DE 49 6605 0101 0075105320

<b>BIC/SWIFT:</b> KARSDE66
<b>BIC/SWIFT:</b> GENODE61KA1
<b>BIC/SWIFT:</b> KARSDE66
<b>BIC/SWIFT:</b> KARSDE66
<b>BIC/SWIFT:</b> KARSDE66

<b>(Konto/account:</b> 226 076 59
<b>(Konto/account:</b> 103 272 10
<b>(Konto/account:</b> 750 041 27
<b>(Konto/account:</b> 750 093 08
<b>(Konto/account:</b> 750 135 40

<b>BLZ:</b> 660 501 01)
<b>BLZ:</b> 661 900 00)
<b>BLZ:</b> 660 501 01)
<b>BLZ:</b> 660 501 01)
<b>BLZ:</b> 660 501 01)